Pulsoksymetr DUO CONTROL

Model: AP 4117

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Infolinia: +48 12 271 33 62

TABLE OF CONTENTS

DESTINY.....	18
INTRODUCTION.....	19
PRODUCT STRUCTURE.....	20
SYMBOLS USED IN THE MANUAL AND ON THE BOX.....	20
PRODUCT FEATURES.....	21
IMPORTANT SAFEGUARDS.....	22
OPERATING YOUR PULSE OXIMETER.....	24
CLEANING AND DISINFECTION.....	27
MAINTENANCE.....	28
TROUBLESHOOTING.....	28
SPECIFICATIONS.....	29
CONTACT.....	30
WARRANTY.....	31

DESTINY

English version

The Pulse Oximeter Sanity, based on all digital technology, is intended for noninvasive spot-check measurement of functional oxygen saturation of arterial hemoglobin (SpO₂). Advanced DSP algorithm* can minimize the influence of motion artifact and improve measurement accuracy of low perfusion*.

The Oximeter can be used to measure human Hemoglobin Saturation and heart rate through finger. The product is suitable for use in family, Oxygen Bar, social medical organizations, physical care in sports and etc.

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*PI (Perfusion Index) is the ratio of the pulsatile blood flow to the non-pulsatile static blood flow in a patient's peripheral tissue, such as finger tip, toe, or ear lobe. Perfusion index is an indication of the pulse strength at the sensor site.

INTRODUCTION

About Sanity®

Sanity® has been present on the Polish market for 20 years as a brand of quality and functional medical products, aiding the return to health, daily health prevention and care for the whole family. Sanity's® product portfolio includes a wide range of medical devices and accessories used to perform basic diagnostic tests, ENT and rehabilitation procedures, as well as treatments supporting daily body care and care for the health and hygiene of adults and children.

Sanity's® products are characterised by their high quality and practical, user-friendly solutions.

Our products comply with the most stringent safety standards, have successfully passed clinical tests and received the corresponding certificates confirming compliance with the EU standards.

Our values

Love is the basis of our humanity. It drives us in life and causes that we become better.


Care is the effect of our feelings and sensitivity. Caring for our close ones guides our actions and changes the way we look at everyday objects.

Help is a sign of love and caring. It includes all, even smallest tasks that we perform for the other person, caring about them in every situation.

PRODUCT STRUCTURE

SYMBOLS USED IN THE MANUAL AND ON THE BOX

	Symbol for "TYPE BF APPLIED PARTS"		Keep dry
IP22	Stopień ochrony zapewniany przez obudowę: IP22		Follow instruction for use
	Important safety instructions		Symbol for „SERIAL NUMBER”
	Symbol for "MANUFACTURER"		Batch code
CE 0197 Symbol for "COMPLIES WITH MDD 93/42/EEC REQUIREMENTS"			
	Symbol for "ENVIRONMENT PROTECTION – Waste electrical products should not be disposed of with household waste. Please recycle where facilities exist. Check with your local authority or retailer for recycling advice".		

PRODUCT FEATURES

- Lightweight for carrying and Easy-To-Use.
- Manually adjust the direction of interface.
- Color OLED display, simultaneous display for testing value and plethysmogram.
- Low Perfusion: 0.2%. (Advanced DSP algorithm can improve measurement accuracy, under the condition of low perfusion.)
- Visual & Sound alarm function. Real-time spot-checks.
- Low Battery voltage indicator.
- Automatically switch off.
- Standard two AAA 1.5V Alkaline Batteries support more than 20 hours continuous work.

Note:

The device can't used to measure the child below 1 year old as the testament result is not guarantee to accurate.

The fingertip pulse oximeter is intended only as an adjunct in patient assessment. It must be used in conjunction with other methods of assessing clinical signs and symptoms.

A function tester cannot be used to assess the accuracy of a pulse oximeter monitor or sensor.

Clinical testing is used to establish the SpO₂ accuracy. The measured arterial hemoglobin saturation value (SpO₂) of the sensor is compared to arterial hemoglobin oxygen(SaO₂) value, determined from blood samples with a laboratory CO-oximeter. The accuracy of the sensors in comparison to the CO-oximeter amples measured over the SpO₂ range of 70 -100%.

PRODUCT FEATURES

Accuracy data is calculated using the root-mean-square (Arms value) for all subjects. Only about two-thirds of PULSE OXIMETER EQUIPMENT measurements can be expected to fall within \pm Arms of the value measured by a CO-Oximeter.

Pulse rate Accuracy shall be used pulse simulator to assess. The measured pulse rate is compared to the preset pulse rate value in simulator. Accuracy data is calculated using the root-mean-square (Arms value) for all subjects.

IMPORTANT SAFEGUARDS

- Do not attempt to service the Pulse Oximeter. Only qualified service personnel should attempt any needed internal servicing.
- Prolonged use or the patient's condition may require changing the sensor site periodically. Change sensor site and check skin integrity, circulatory status and correct alignment at least every 2 hours.
- SpO₂ measurements may be adversely affected in the presence of high ambient light. Shield the sensor area (with a surgical towel, or direct sunlight, for example) if necessary.
- The following reason will cause interference to the testing accuracy of the pulse oximeter.
 - a) High-frequency electrosurgical equipment.
 - b) Placement of a sensor on an extremity with a blood pressure cuff arterial catheter, or intravascular line
 - c) The patient has hypotension severe vasoconstriction severe anemia or hypothermia.

IMPORTANT SAFEGUARDS

- d) The patient is in cardiac arrest or is in shock.
- e) The patient is in cardiac arrest or is in shock.
- f) Fingernail polish or false fingernails may cause inaccurate SpO₂ readings.

Caution!

EXPLOSION HAZARD – Do not use the Pulse Oximeter in a flammable atmosphere where concentrations of flammable anesthetics or other materials may occur.

Do not throw batteries in fire as this may causes them to explode.

Do not attempt to recharge normal dry-cell batteries, they may leak. And may cause a fire or even explode.

Do not use the pulse oximeter in an MRI or CT environment.

Do not modify this equipment without authorization of the manufacturer.

If this equipment is modified, appropriate inspection and testing must be conducted to ensure continued safe use of equipment.

Keep the operating environment free of dust, vibrations, corrosive, or flammable materials, and extremes of temperature and humidity.

Do not operate the unit if it is damp or wet because of condensation or spills. Avoid using the equipment immediately after moving it from a cold environment to a warm, humid location.

Never use sharp or pointed objects to operate the front-panel switches.

The batteries must be taken out from the battery compartment if the device will not be used for a long time.

IMPORTANT SAFEGUARDS

The device shall only be used if the battery cover is closed.

The batteries must be properly disposed according to local regulation after their use.

The device should be kept away from children, pets and pests to avoid swallowing.

OPERATING YOUR PULSE OXIMETER

Parameter setting

Press the direction button for one second, the oximeter will enter into parameter setting. There are two submenus for choice:

When the * signal is shown on the Alarm Setup, press the direction button for 1s and enter into the Sound Setting menu (fig.1), press the direction button for 1s to turn on/off for the alarm and beep.

When the * signal is shown on the Sounds Setup, press the direction button for 1s and enter into the Alarm Setting menu (fig.2), you can press the direction button in turn to select the item. And press the direction button for 1s to change the data you need.

OPERATING YOUR PULSE OXIMETER

On the Sounds Setup menu (fig.1), when the * signal shown on front of the "+/-" sign "+" on the right side can be changed to "-", by pressing the direction button for 1s. When "+" shows on the right side, press the direction button for 1s can increase the SPO2 and PR's highest and lowest alarm value. When "-" shows on the right side, press the direction button for

1s can reduce the SPO2 and PR's highest and lowest alarm value.

On the Alarm Setup menu (fig.2). When the * signal shown on "Restore", press the direction button for 1s, the right side can be changed to "OK", which cause the device to restore factory setting data.

Fig.1

Fig.2

Note:

The alarm have 1 second delay after the incorrect result being detected. The customer can preset the alarm value to the 98 or 99 to check whether it is normal for alarm setting.

Install battery

- Turn the device upside down and open the battery cover.
- Installing two AAA batteries into battery cassette in correct polarities and cover it

WARNING! Do not attempt to recharge normal alkaline batteries, they may leak and may cause a fire or even explode.

OPERATING YOUR PULSE OXIMETER


Turn the Pulse Oximeter on

Put one of fingers into rubber hole of the oximeter (it is best to put the finger thoroughly) with nail surface upward, then releasing the clamp.

Press power button to turn the Pulse Oximeter on. The oximeter will automatically be powered off when no finger in the device for longer than 16 seconds.

Display Description of OLED

The display interface of "OLED" can rotate four directions with six different display modes after pressing the power button for less than 0.5s. It is shown as below:

Note:

1. When battery power is at lowest level, the battery capacity indicates symbol of  in OLED, remind users of replacement of battery.
2. The pletismograph can be regarded as correct if the wave is fluctuated regularly.

CLEANING AND DISINFECTION

Cleaning

Cleaning the machine once every day. Switch off the power and take out the batteries before cleaning. Cleaning exterior surface (OLED display screen included) of the unit with a dry and soft cloth. Use 75% density of medical alcohol to clean the surface and use dry fabric with little alcohol to avoid alcohol permeates into the device.

Disinfection

Disinfect the machine after using by the patient if multiple patients use the machine in the hospital.

Use 75% density of medical alcohol to clean the surface that contacting with the patient.

Caution!

Don't use strong solvent. For example, acetone.

Never use an abrasive such as steel wool or metal polish.

Do not allow any liquid into the product, and do not immerse any parts of the device into any liquids.

Avoid pouring liquids on the device while cleaning.

Don't remain any cleaning solution on the surface of the device.

MAINTENANCE

- Replace the batteries timely when battery indication is low. Clean surface of the Pulse Oximeter before it is used in diagnosis for patients.
- Remove the batteries inside the battery cassette if the oximeter will not be operated for a long time.
- It is better to preserve the product in a place where ambient temperature is -10-40C and humidity is 10%-80%.
- Regular inspection to make sure that no obvious damage existed to affect the safety and performance of device.
- No flammable substance, otoprot or lower temperature and humidity existed in operation conditions.

TROUBLESHOOTING

PROBLEMS	POSSIBLE REASON	RESOLUTIONS
Oxyhemoglobin or heart rate can not be shown normally	<ol style="list-style-type: none">1. Finger is not plugged correctly.2. Patient's perfusion is too low to be measured.	<ol style="list-style-type: none">1. Retry by plugging the finger2. Try some more times, if you can make sure about no problem existing in the product, Please go to a hospital timely for exact diagnosis
Oxyhemoglobin or heart rate is shown unstably	<ol style="list-style-type: none">1. Finger might not be plugged deep enough2. Finger is trembling or patient's body is in movement status	<ol style="list-style-type: none">1. Retry by plugging the finger2. Try not to move, Let the patient keep calm.

TROUBLESHOOTING

Oxyhemoglobin or heart rate is abnormal and cause alarm	<ol style="list-style-type: none">1. Finger is not plugged correctly.2. Patient's SPO2&PR is abnormal.	<ol style="list-style-type: none">1. Retry by plugging the finger2. go to the hospital for further examination.
The oximeter cannot be powered on	<ol style="list-style-type: none">1. Power of batteries might be inadequate or not be there at all2. Batteries might be installed incorrectly3. The Oximeter might be damaged	<ol style="list-style-type: none">1. Please replace batteries2. Please reinstall the batteries3. Please contact with local customer service center
The screen are suddenly off	<ol style="list-style-type: none">1. The product is automatically powered off when no signal is detected longer than 8 seconds2. Power quantity of the batteries is exhausted.	<ol style="list-style-type: none">1. Normal2. Replace the batteries

SPECIFICATIONS

Physical Characteristics	Machine: Dimensions -63 mm (L) x 36mm (W) x 34mm (D) Weight -approx: 35 g Outer box: Dimensions: 87 mm (L) x 68 mm (W) x 40 mm (D), Weight: 85g
Anti-electric Shock Type:	Internally powered equipment
Anti-electric Shock Degree	Type BF equipment
EMC	Type A class II
Mode of operation:	Continuous Operation.
Enclosure Degree of ingress protection:	IP22

SPECIFICATIONS

Power	Internal: 2xAAA 1.5v alkaline battery Power Consumption: Smaller than 30mA (Normal)
Environmental	Operating Temperature: 10°C do 40°C Storage Temperature: -20°C do 70°C Relative Humidity: 15% do 90% (non-condensing) Air Pressure: 86Kpa-106Kpa Operating Attitude: 0-2000 m
Alarm default value:	Hemoglobin saturation: Upper limit: 100/ bottom limit:90 Pulse rate: Upper limit: 130 /bottom limit: 50
Electronics Parameters:	Hemoglobin saturation display: 35-100% Pulse rate Display: 30-250 BPM Resolution: – Hemoglobin Saturation 1% – Pulse rate 1 BPM
Measure Accuracy:	Hemoglobin Saturation: $\pm 2\%$ (70%-100%) Unspecified (<70%) Pulse rate: ± 2 BPM
Probe LED Specification	RED: Wave Length: 660 ± 2 nm Radiant Power: 1,8 mW Infra RED: Wave Length: 905 ± 2 nm Radiant Power: 2,0 mW

CONTACT

For more information about our products, please visit www.sanity-products.com

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